



Pen Needles

Trusted performance at a competitive price

With savings up to 60% versus other leading national brands¹, Unifine® delivers trusted performance at a competitive price. As one of the most trusted pen needle brands in the U.S., we put the patient first with innovative features that help enhance quality of life and drive adherence. We make patient access a priority by investing in programs that improve accessibility and affordability. Unifine® pen needles are covered through most insurance plans² and are available in a broad range of sizes to meet unique needs and preferences.



 **Unifine® SafeControl®**



 **Unifine® Pentips® Plus**



 **Unifine® Pentips®**

OWEN MUMFORD

1755 West Oak Commons Court | Marietta, GA 30062
1-800-421-6936 | owenmumford.com/us

Take CONTROL

Other safety pen needles may **prematurely activate** before the dose is delivered, Unifine® SafeControl® is designed to provide the right balance of **safety** and **control**, from injection to retraction.

- Dual needle protection to help prevent needlestick injuries
- Enhanced user control preventing premature activation³
- Visual confirmation of dose-delivery³
- Ideal for use in long-term care
- Covered by most Medicare Part D plans



Before activation

Visible needle
utilizes the same
technique as
conventional pen
needles



**Innovative
push tab**
for one-handed
control of the safety
mechanism



Green cap
for sterility and
protection from
the needle prior
to injection




After activation


Safety guard
for covering the needle
after injection to help
prevent needlestick
injuries



**Dual safety
indicator**
for visual and audible
confirmation



Needle covered
at both ends to help
avoid needlestick
injuries



Visual confirmation of dose-delivery

Intuitively designed to provide enhanced user control.³



Universal Fit

Compatible with most injection pens⁴



Made with precision in the UK

Unifine® Pentips® Plus

Pen Needle with a Built-In Remover

Preferred by patients^{4,5},
Unifine® Pentips® Plus is
designed to make pen
needle removal **easier**,
safer, and more
convenient⁵.

In an eight week observational behavioral study, patients felt Unifine Pentips Plus was easy to use, offered added convenience and a safer disposal option^{5,6}

61% of users preferred
Unifine Pentips Plus^{4,5}

86% of users found Unifine Pentips
Plus easy to change^{5,7}


59% of users felt Unifine Pentips Plus was
more convenient when on-the-go^{5,7}

61% increase in patient adherence
versus standard pen needle^{5,9}


OptiFlow®
TECHNOLOGY

Designed for improved drug flow and
low delivery force⁸

With OptiFlow®



Without OptiFlow®



Unifine® Pentips®

Pen Needle

Tried, true and **trusted**.

Unifine® Pentips® deliver the same European-manufactured quality and exacting standards as Unifine® Pentips® Plus. Available in a broad range of sizes.

 **DiamondPoint**
TECHNOLOGY™

Sharpened to perfection for
lower penetration force⁸

 **LowGlideForce**

Electropolished and treated with advanced
silicone for exceptional comfort

UniversalFit 







Compatible with most injection pens⁴

PRODUCT RANGE **Flexible options** Available in a broad range of sizes and packaging options (100ct & 30ct).

Unifine SafeControl®
Safety Pen Needle with Dual-Protection

					
4mm x 32G	5mm x 31G	5mm x 30G	6mm x 31G	8mm x 31G	8mm x 30G
100 COUNT Item # AN 7940 UPC 384707940013 RC* 08470-7940-01	100 COUNT Item # AN 7950 UPC 384707950012 RC* 08470-7950-01	100 COUNT Item # AN 7955 UPC 384707955017 RC* 08470-7955-01	100 COUNT Item # AN 7990 UPC 384707990018 RC* 08470-7990-01	100 COUNT Item # AN 7930 UPC 384707930014 RC* 08470-7930-01	100 COUNT Item # AN 7935 UPC 384707935019 RC* 08470-7935-01

Unifine Pentips® Plus
Pen Needle with a Built-In Remover

					
4mm x 33G	4mm x 32G	5mm x 31G	6mm x 31G	8mm x 31G	12mm x 29G
100 COUNT Item # AN 3860 UPC 384703860018 RC* 08470-3860-01	100 COUNT Item # AN 3840 UPC 384703840010 RC* 08470-3840-01	100 COUNT Item # AN 3850 UPC 384703850019 RC* 08470-3850-01	100 COUNT Item # AN 3890 UPC 384703890015 RC* 08470-3890-01	100 COUNT Item # AN 3830 UPC 384703830011 RC* 08470-3830-01	100 COUNT Item # AN 3829 UPC 384703829015 RC* 08470-3829-01
10% Thinner than 32G ¹⁰	30 COUNT Item # AN 1840 UPC 384701840012 RC* 08470-1840-01	30 COUNT Item # AN 1850 UPC 384701850011 RC* 08470-1850-01	30 COUNT Item # AN 1890 UPC 384701890017 RC* 08470-1890-01	30 COUNT Item # AN 1830 UPC 384701830013 RC* 08470-1830-01	

Unifine Pentips®
Pen Needle

						
4mm x 33G	4mm x 32G	5mm x 31G	6mm x 32G	6mm x 31G	8mm x 31G	12mm x 29G
100 COUNT Item # AN 3560 UPC 384703560017 RC* 08470-3560-01	100 COUNT Item # AN 3540 UPC 384703540019 RC* 08470-3540-01	100 COUNT Item # AN 3550 UPC 384703550018 RC* 08470-3550-01	100 COUNT Item # AN 3595 UPC 384703595019 RC* 08470-3595-01	100 COUNT Item # AN 3590 UPC 384703590014 RC* 08470-3590-01	100 COUNT Item # AN 3530 UPC 384703530010 RC* 08470-3530-01	100 COUNT Item # AN 3529 UPC 384703529014 RC* 08470-3529-01
10% Thinner than 32G ¹⁰	30 COUNT Item # AN 1140 UPC 384701140013 RC* 08470-1140-01	30 COUNT Item # AN 1150 UPC 384701150012 RC* 08470-1150-01		30 COUNT Item # AN 1190 UPC 384701190018 RC* 08470-1190-01	30 COUNT Item # AN 1130 UPC 384701130014 RC* 08470-1130-01	

*Reimbursement Code. 1. Comparison of acquisition cost of Unifine Pentips versus other leading national brands. 2. Independent and chain pharmacy co-pay adjudication (April 2016). 3. Project Saturn B (2017) Face-to-face interviews and product evaluations commissioned with an independent market research agency. Data on file. 4. Data on file. 5. HRW (2014) Impact of Unifine Pentips Plus on pen needle changing behaviour amongst people with diabetes medicating with injectable formats. 6. The used pen needle is secured in the interim container and must be placed in a sharps container at the earliest opportunity. 7. Versus standard pen needle. 8. Melab GmbH Testing 2017. 9. Adherence is defined as the changing of pen needle after every injection. 10. Comparison of 33G to 32G max outer needle diameter in accordance with ISO 9626:2016 standards.

UBR/OMI/0824/1/US

OWEN MUMFORD

1755 West Oak Commons Court | Marietta, GA 30062
1-800-421-6936 | owenmumford.com/us